

New products  
In worldwide

# 「Vital Signature」 Sphygmomanometer

- Korotokoff method of the world Medical Science Association standard
- The process of the blood pressure value measurement is displayed.
- The shock state sphygmomanometry is possible.
- The blood pressure of the newborn baby can be measured.
- The sphygmomanometry when moving is possible.

\*Patent pending

## Existing sphygmomanometer

- Mercury type:  
The only product that Medical Science Association admits including WHO. The computerization of the measurement data cannot be done. The individual variation by the measurer goes out.
- Cuff oscillometric methods:  
It is not used for the diagnosis and treatment



## Vital Signature Sphygmomanometer


Electronic measurement type that adopted korotokoff method that Medical Science Association in the world admits. It is possible to display the proof reason for the measurement data in real time.



Sphygmomanometry



Main body and touch panel PC




Measuring corrugate display

Automation of measurement



Real-time display of proof data



Monitor indication  
An interlocking movement with the electronic chart


Print out




FDA in Japan

Synchronization with 200 points of measurement medical treatment fee of 24 hours.(in Japan)

Vital Signature  
Sphygmomanometer series



The shock state can be measured.



The blood pressure  
Of the newborn baby can be measured.



Measurement when moving



Healthcare network